

wellsystem

INFOCARD HYDROWAVE BED

THE FIRST DESIGNER FURNITURE PIECE FOR WELL-BEING AND RELAXATION





HydroWave Bed

Soothing dry water massage meets gentle infrared warmth - in the proven quality of **well**system. Elegantly designed, fully customizable, and ideal for exclusive interiors.

Just 15 minutes on the **well**system HydroWave Bed noticeably enhances well-being and provides pure relaxation. Simply lie back and enjoy.

| | | |
|-----------------------|---|---|
| OPERATION | Operation | Quick-start operation and smartphone app (iPhone & Android) |
| | Professional set-up for operator settings | No |
| | Massage programs can be stored | Massage routines can be saved via the app |
| MASSAGE | No credit system | Start directly on the device or through the app |
| | Massage zones | Full body, Upper body, Lower body |
| | Massage techniques | Stroking massage, Mirror massage, Circular massage, Zigzag massage |
| | Pre-installed massage programs | Pure Relaxation, Energy Boost, Healthy Back, Fit Legs |
| | Pump performance | 6,5 bar |
| | Selectable massage intensity | Multi-level |
| | Adjustable massage speed | No |
| | Adjustable body length | Yes |
| | Size of the lying surface | 2.080 x 1.190 mm |
| | Water cooling (optional) | No |
| COOLING | Air cooling | No |
| DESIGN | Device color | Various color combinations available: Upholstery: anthracite or beige Decorative panel: white, black, or walnut Base: white or black |
| TECHNICAL INFORMATION | Device dimensions (L x W x H) | 2.580 x 1.190 x 770 mm |
| | Weight (empty unit) | 168 kg |
| | Maximum user weight | 150 kg |
| | Maximum session duration | 60 minutes 150 minutes per day |
| | Temperature | 25–50 °C, adjustable infrared heat temperature |
| | Water volume | approx. 304 liters |
| | Maximum power consumption | 2.900 watt |
| | Art. No | Article number depends on the selected color combination |